



UNIPLAST

Different By Design



Plastic Lab Ware Catalogue

2025-26



Celebrating our Past
Inspiring the Future



Founded in 1969 in New Delhi, UNIPLAST is a trusted manufacturer of high-quality laboratory and scientific plasticware, proudly serving the scientific, pharmaceutical, and cosmetic industries for over 55 years. Backed by ISO 9001:2008 certification and a team of qualified engineers, we blend deep manufacturing expertise with relentless quality control.

With advanced in-house capabilities for injection and blow molding, we craft a wide array of labware using premium-grade polymers. From microcentrifuge tubes and inoculation loops to amber bottles, desiccators, and specimen jars, every product reflects our commitment to precision, durability, and elegant design.

Our infrastructure includes clean room environments, modern tooling systems, and a robust knowledge base—enabling us to deliver custom solutions in varied sizes, designs, and colors. UNIPLAST products are distinguished by consistent branding, ergonomic shapes, and crystal-clear calibrations that support scientific accuracy and visual clarity.

We don't just manufacture — we empower innovation. This catalog reflects our diverse product line and commitment to continuous improvement, keeping pace with emerging research and end-user needs.

Shiksha Tyagi

Director



Marketing Office: BP 174, Shalimar Bagh (West), Delhi 110088, India

Production Facility: D 221, Sector 3, D.S.I.I.D.C., Bawana Ind. Area, New Delhi - 110039, India

Tele: +91-9971700758, +91-9971700756, +91-8700221715, +91-7827065011, +91-11-45520870

Email: UNIPLAST@YMAIL.COM , UNIPLAST@UNIPLAST.US

Website: Uniplast.org.in

Index

Bottles

Dropper Bottles T. K	01
Dropper Bottles with Tip Cap	01
Narrow Mouth Bottles	04
Wide Mouth Bottles	05
Sterile Wide Mouth Bottle	06
Print Graduation - Wide Mouth Bottle	06
Print Graduation - Narrow Mouth Bottle	06
Conical Flask Printed / Engraved	07
Volumetric Flask	07
Rectangular Wide Mouth Bottle	07
Large Wide Mouth Bottle with Handle	08
Large Narrow Mouth Bottle	08
Chemical Wide Mouth Bottles	10
Oblong Jar with Handle	10
Squeeze (Gel) Bottle Dispenser	10
Pharmaceutical Powder & Tablet Container	10

Wash Bottles

Wash Bottles - Regular & With Dispenser	01
Wash Bottles Wide Mouth Self-Venting	02
Wash Bottles Narrow Mouth Self-Venting	02
Wash Bottles Wide Mouth Autoclavable	02
Wash Bottles Narrow Mouth Autoclavable	02
Self-Venting Labelled Wide Wash Bottle	06
Self-Venting Labelled Narrow Wash Bottle	06

Carboys

Aspirator Bottles with Stopcock	09
Carboys	09
Carboys Wide Mouth	09
Carboys with Spigot	09

Volumetrics

Measuring Cylinders - Hexa and Round Base	11
Measuring Jug	11
Measuring Beaker	11
Measuring Jugs - Printed	12
Measuring Beaker - Printed	12
Rain Gauze Cylinder	12
Rain Gauze Set	12

Desiccators

Desiccator (Vacuum)	15
Desiccator (Plain)	15
Spare Silicon 'O'-Ring for Desiccator	15

Fluid Handling

Connectors "T", "Y", "L"	30
--------------------------	----

Funnels

Funnel Holder	14
Separatory Funnel Holder	14
Separatory Funnel	14
Analytical Funnel	14

Hospital Consumables

Kidney Tray	18
Bed Pan	18
Urine Pot (Male and Female)	18
Sputum Mug with Lid	19
Sputum Specimen Container	19
Stool Specimen Container	19
Sample Container for Urine	19
Laboratory Jars with Screw Caps	20

Micro Biology Apparatus

L Shaped Spreader	16
Soft Loops	16
Stirring Rods	16
Spatula	17
Measuring Scoop	17
Swab Stick	20

Liquid Handling

Micro Pipettes Fixed Volume	25
Micro Pipettes Variable Volume	25
Pipettor Stand - Z Type	25
Micro Pipette Stand with Drawer	25
Bottle Top Dispenser	26
Micropipette Tips and Filter Tips	26
Pipette Tip Box	26
Pipette Pump	26
Pasture Pipette	27
ESR Pipettes (Disposable)	27
ESR Pipette Stand	27

Tissue Culture

Petri Dish Autoclavable	15
Petri Dish Disposable - Sterilized	15
Rack for Petri Dish	16
Specimen Jar / Gas Jar	20
Wide Mouth Specimen Tube	20

Education and Atomic Models

Animal Cage	28
Plantation Pots	28
Kipps Apparatus	28
Over Flow Jar	28
Paschal Law Apparatus	29
Pneumatic Trough	29
Simple Cell Pot	29
Leclanché Cell Pot	29
Siphon Automatic Apparatus	30
Crystal Model Set	31
Atomic Model Set	31

Stands and Racks

Retort Stand	13
Burette Clamp (Single/Double)	13
Burette Holder	13
Pipette Stand (Vertical)	24
Pipette Rack (Horizontal)	24
Stand for MCT, RIA & Storage Vial	23
Test Tube Stand (3 Tier)	23
Test Tube Stand 16 x25	23
Round Test Tube Stand	24
Draining Rack	24
Conical Centrifuge Tube Rack	22
Laboratory Tray	27

Slide Accessories

Slide Box Detachable and Hinged	17
Slide Box Hinged	17
Slide Mailer	17
Coplin Jar	18

Vials

Specimen Backman & Hitachi Cup	21
Ria Vial	21
Screw Top Micro Centrifuge tube	21
Micro Centrifuge Tube	21
Centrifugal Tubes Conical Bottom	22
Centrifugal Tubes Self-Standing	22
Storage Vials (self-standing)	22
Cryo Vial	22
Cryo Box	23

01 Dropping Bottles

Material: Low Density Polyethylene
Colour: Natural White, Translucent

HSN Code: 39233090

Dropping bottles are used for repeated, consistent, accurate dispensing of reagents, one drop at a time or in a stream.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101201	30 ml.	12 Pcs.	14.00	168.00
101202	60 ml.	12 Pcs.	17.00	204.00
101203	125 ml.	12 Pcs	21.00	252.00



02 Dropping Bottle with Tip Cap

Material: Low Density Polyethylene
Colour: Natural White, Translucent

HSN Code: 39233090

Dropper bottles are transparent with an extra-long, fine dropping tip for accurate dispensing.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101305	30 ml.	12 Pcs.	15.50	186.00
101306	60 ml.	12 Pcs.	18.00	216.00
101307	125 ml.	12 Pcs	22.00	264.00
101308	250 ml.	12 Pcs	34.50	414.00



03 Wash Bottles

Material: Low Density Polyethylene with PP Closure
Colour: Natural White, Translucent (Closures available in multi colour)

HSN Code: 39233090

These bottles are translucent and flexible yet unbreakable. Solutions can be dispensed easily by lightly squeezing the bottle, used to rinse various laboratory wares. Wash bottles are sealed with a screw-top lid. The flexible polyethylene delivery tube can be directed as appropriate.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101401	125 ml.	12 Pcs.	25.00	300.00
101402	250 ml.	12 Pcs.	33.50	402.00
101403	500 ml.	12 Pcs.	43.50	522.00
101404	1000 ml.	06 Pcs.	75.00	450.00



04 Wash Bottles Side Dispenser

Material: Low Density Polyethylene with PP Closure
Colour: Natural White, Translucent (Closures available in multi colour)

HSN Code: 39233090

These translucent bottles ensure a leak-proof design and feature a wide mouth for easy filling. They are flexible yet unbreakable. The standout feature is the built-in delivery tube, making it a one-piece wash bottle. Ideal for general-purpose solutions such as water, soap, dilute acids, and dilute bases.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101501	125 ml.	12 Pcs.	36.50	438.00
101502	250 ml.	12 Pcs.	52.00	624.00
101503	500 ml.	12 Pcs.	67.50	810.00





1000 and 500 ml size – Wide Mouth - LDPE

05 Wide Mouth Wash Bottle (LDPE with PP closure and Draw Tube)

Material: Low Density Polyethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

These self-venting bottles moulded in LDPE are Soft and easy to squeeze. Plain without label can be used universally as per chemical needs. Feature graduated markings for easy measurement. The closure and stem are moulded in one piece, ensuring a secure, leak-proof seal. The tip can be cut back to increase flow.

Code	Capacity	Description	Packing	Price / Piece (₹)
101705	500 ml.	Plain Bottle with Plain Closure	06 Pcs.	210.00
101705C	500 ml.	Plain Bottle with Coloured Closure	06 Pcs.	230.00
101805	1000 ml.	Plain Bottle with Plain Closure	06 Pcs.	270.00
101805C	1000 ml.	Plain Bottle with Coloured Closure	06 Pcs.	290.00
SV1018	Spare cap	WM Screw/ Stem/Draw Tube	01 Pcs.	210.00



250 and 500 ml size – Narrow Mouth - LDPE

06 Narrow Mouth Wash Bottle (LDPE with PP closure and Draw Tube)

Material: Low Density Polyethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

These self-venting Bottles are “Plain- without label” and can be used universally as per the chemical needs of the lab. The Closures are also available in different colours The closure and stem are moulded in one piece for a secure leak proof seal.

Code	Capacity	Description	Packing	Price / Piece (₹)
101605	250 ml.	Plain Bottle with Plain Closure	06 Pcs.	180.00
101605C	250 ml.	Plain Bottle with Coloured Closure	06 Pcs.	200.00
101905	500 ml.	Plain Bottle with Plain Closure	06 Pcs.	200.00
101905C	500 ml.	Plain Bottle with Coloured Closure	06 Pcs.	220.00
SV1016	Spare cap	NM Screw/ Stem/Draw Tube	01 Pcs.	190.00



500 and 1000 ml size - (Autoclavable) Polypropylene

07 Wide Mouth Autoclavable Wash Bottle - PP

Material: Polypropylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

Bottle is moulded in Polypropylene Specifically designed to be Soft, squeezable and suitable for repeated autoclaving for pharma, micro lab. The closure and stem are moulded in one piece, ensuring a secure, leak-proof seal The dispensing tip can be cut back to increase flow. Ideal for general-purpose solutions like saline and buffers that may require autoclaving.

Code	Capacity	Description	Packing	Price / Piece (₹)
101707	500 ml.	Plain Bottle with Plain Closure	06 Pcs.	210.00
101707C	500 ml.	Plain Bottle with Coloured Closure	06 Pcs.	230.00
101807	1000 ml.	Plain Bottle with Plain Closure	06 Pcs.	270.00
101807C	1000 ml.	Plain Bottle with Coloured Closure	06 Pcs.	290.00
SV1018	Spare cap	WM Screw/ Stem/Draw Tube	01 Pcs.	210.00



250 and 500 ml size - (Autoclavable) Polypropylene

08 Narrow Mouth Autoclavable Wash Bottle - PP

Material: Polypropylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

Soft and squeezable. Specifically designed for repeated autoclaving for pharma, micro lab. The closure and stem are moulded in one piece, ensuring a secure, leak-proof seal. The dispensing tip can be cut back to increase flow. These self-venting Bottles are “Plain- without label” and can be used universally as per the chemical needs of the lab. The Closures are also available in different colours.

Code	Capacity	Description	Packing	Price / Piece (₹)
101607	250 ml.	Plain Bottle with Plain Closure	06 Pcs.	180.00
101607C	250 ml.	Plain Bottle with Coloured Closure	06 Pcs.	200.00
101907	500 ml	Plain Bottle with Plain Closure	06 Pcs.	200.00
101907C	500 ml.	Plain Bottle with Coloured Closure	06 Pcs.	220.00
SV1016	Spare cap	NM Screw/ Stem/Draw Tube	01 Pcs.	190.00

09 Self-Venting Labelled Wide Wash Bottle - LDPE

Material: Low Density Polyethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

Each self-venting bottle is designed with a one-piece moulded closure, stem & draw tube, ensuring a secure, leak-proof seal. A colour bar matching the closure helps prevent cross-contamination. These non-metallic, self-venting wash bottles are clearly labelled with chemical name, formula, & CAS no.

09 A 500 ml Size – Wide Mouth (Labelled)- LDPE with PP Closure & Draw Tube

Code	Capacity	Chemical Name	Packing	Price / Piece (₹)
101710	500 ml.	Acetone	06 Pcs.	175.00
101711	500 ml.	Distilled Water	06 Pcs.	175.00
101712	500 ml.	Iso Propanol	06 Pcs.	175.00
101713	500 ml.	Ethanol	06 Pcs.	175.00
101714	500 ml.	Methanol	06 Pcs.	175.00
101715	500 ml.	Sodium Hypo Chlorite	06 Pcs.	175.00
101716	500 ml.	Mix Pack (6 above Chemicals)	06 Pcs.	200.00
101717	500 ml.	Ethyl Acetate	06 Pcs.	175.00
101718	500 ml.	Toluene	06 Pcs.	175.00
101719	500 ml.	Methyl Ethyl Ketone	06 Pcs.	175.00
101720	500 ml.	Xylene	06 Pcs.	175.00



09 A

500 ml Size – LDPE (Labelled)

09 B 1000 ml Size – Wide Mouth (Labelled)- LDPE with PP Closure & Draw Tube

Code	Capacity	Chemical Name	Packing	Price / Piece (₹)
101810	1000 ml.	Acetone	06 Pcs.	310.00
101811	1000 ml.	Distilled Water	06 Pcs.	310.00
101812	1000 ml.	Iso Propanol	06 Pcs.	310.00
101813	1000 ml.	Ethanol	06 Pcs.	310.00
101814	1000 ml.	Methanol	06 Pcs.	310.00
101815	1000 ml.	Sodium Hypo Chlorite	06 Pcs.	310.00
101816	1000 ml.	Mix Pack (6 above Chemicals)	06 Pcs.	340.00
101817	1000 ml.	Ethyl Acetate	06 Pcs.	310.00
101818	1000 ml.	Toluene	06 Pcs.	310.00
101819	1000 ml.	Methyl Ethyl Ketone	06 Pcs.	310.00
101820	1000 ml.	Xylene	06 Pcs.	310.00
SV1018	Spare cap	WM Screw/ Stem/Draw Tube	01 Pcs.	210.00



09 B

1000 ml Size – LDPE (Labelled)

10 Self-Venting Labelled Narrow Wash Bottle - LDPE

Material: Low Density Polyethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent (Closures available in multi colour)

The closure and stem are moulded in one piece for a secure leak proof sea, stem and draw tube. Each bottle has a colour bar which corresponds to its coloured closure to prevent cross contamination.

10 A 250 ml Size – Narrow Mouth (Labelled)- LDPE with PP Closure & Draw Tube

Code	Capacity	Chemical Name	Packing	Price/ Piece (₹)
101610	250 ml.	Acetone	06 Pcs.	220.00
101611	250 ml.	Distilled Water	06 Pcs.	220.00
101612	250 ml.	Iso Propanol	06 Pcs.	220.00
101613	250 ml.	Ethanol	06 Pcs.	220.00
101614	250 ml.	Methanol	06 Pcs.	220.00
101615	250 ml.	Sodium Hypo Chlorite	06 Pcs.	220.00
101616	250 ml.	Mix Pack (6 above Chemicals)	06 Pcs.	250.00
101617	250 ml.	Ethyl Acetate	06 Pcs.	220.00
101618	250 ml.	Toluene	06 Pcs.	220.00
101619	250 ml.	Methyl Ethyl Ketone	06 Pcs.	220.00
101620	250 ml.	Xylene	06 Pcs.	220.00



10 A

250 ml Size – LDPE (Labelled)

10 B 500 ml Size – Narrow Mouth (Labelled)- LDPE with PP Closure & Draw Tube

Code	Capacity	Chemical Name	Packing	Price/ Piece (₹)
101910	500 ml.	Acetone	06 Pcs.	240.00
101911	500 ml.	Distilled Water	06 Pcs.	240.00
101912	500 ml.	Iso Propanol	06 Pcs.	240.00
101913	500 ml.	Ethanol	06 Pcs.	240.00
101914	500 ml.	Methanol	06 Pcs.	240.00
101915	500 ml.	Sodium Hypo Chlorite	06 Pcs.	240.00
101916	500 ml.	Mix Pack (6 above Chemicals)	06 Pcs.	270.00
101917	500 ml.	Ethyl Acetate	06 Pcs.	240.00
101918	500 ml.	Toluene	06 Pcs.	240.00
101919	500 ml.	Methyl Ethyl Ketone	06 Pcs.	240.00
101920	500 ml.	Xylene	06 Pcs.	240.00
SV1016	Spare cap	NM Screw/ Stem/Draw Tube	01 Pcs.	190.00



10 B

500 ml Size – LDPE (Labelled)


11

11 Narrow Mouth Reagent Bottle – PP

Material: Polypropylene (In Conformance with US FDA 21 CF)

HSN Code: 39233090

Colour: Natural, Transparent

These bottles are leakproof, rigid, translucent & has excellent contact clarity. Ideal for storing & transporting chemicals. Can be autoclaved (Set the closure on top of the container without engaging the threads before autoclaving)

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
102100	4 mL	120 Pcs.	5.75	690.00
102101	8 mL	120 Pcs.	7.00	750.00
102102	15 mL	120 Pcs.	7.50	870.00
102103	30 mL	120 Pcs.	9.75	1230.00
102104	60 mL	120 Pcs.	12.75	1590.00
102105	125 mL	96 Pcs.	16.50	1632.00
102106	250 mL	48 Pcs.	28.00	1344.00
102107	500 mL	24 Pcs.	35.50	828.00
102108	1000 mL	18 Pcs.	69.50	1251.00


12

12 Narrow Mouth Reagent Bottle – HDPE

Material: High Density Poly Ethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent

Low in trace elements. Leak-proof. Rigid with excellent chemical resistance to most acids, bases and alcohols. Ideal for storing, packaging, and shipping liquids. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
102200	4 mL	120 Pcs.	5.75	690.00
102201	8 mL	120 Pcs.	6.00	720.00
102202	15 mL	120 Pcs.	7.00	840.00
102203	30 mL	120 Pcs.	10.00	1200.00
102204	60 mL	120 Pcs.	13.00	1560.00
102205	125 mL	96 Pcs.	17.00	1632.00
102206	250 mL	48 Pcs.	28.00	1344.00
102207	500 mL	24 Pcs.	34.50	828.00
102208	1000 mL	18 Pcs.	69.50	1251.00


13

13 Narrow Mouth Reagent Bottle – Amber

Material: High Density Poly Ethylene with PP Closure

HSN Code: 39233090

Colour: Amber, Opaque

Leak-proof & designed to protect light-sensitive liquids. Rigid with excellent chemical resistance to most acids, bases & alcohols. Ideal for storing, packaging & shipping. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
102300	4 mL	120 Pcs.	5.75	690.00
102301	8 mL	120 Pcs.	6.50	780.00
102302	15 mL	120 Pcs.	7.75	930.00
102303	30 mL	120 Pcs.	11.75	1410.00
102304	60 mL	120 Pcs.	14.00	1680.00
102305	125 mL	96 Pcs.	17.50	1680.00
102306	250 mL	48 Pcs.	29.00	1392.00
102307	500 mL	24 Pcs.	36.50	876.00
102308	1000 mL	18 Pcs.	74.00	1332.00


14

14 Narrow Mouth Reagent Bottle - LDPE

Material: Low Density Poly Ethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent

These translucent bottles are leakproof having. Chemical Resistance, excellent flexibility yet can resist impact. Contents are easily visible. These Bottles are not autoclavable. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
102401	8 mL	120 Pcs.	6.50	780.00
102402	15 mL	120 Pcs.	7.25	870.00
102403	30 mL	120 Pcs.	10.50	1260.00
102404	60 mL	120 Pcs.	13.00	1500.00
102405	125 mL	96 Pcs.	17.50	1584.00
102406	250 mL	48 Pcs.	28.00	1296.00
102407	500 mL	24 Pcs.	34.00	816.00
102408	1000 mL	18 Pcs.	69.00	1242.00

13 Wide Mouth Reagent Bottle – PP

Material: Polypropylene (In Conformance with US FDA 21 CF)

HSN Code:39233090

Colour: Natural, Transparent

These translucent round bottles are leakproof having excellent flexibility yet can resist impact. Contents are easily visible. These Bottles are not autoclavable.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
102603	30 ml.	120 Pcs.	11.50	1350.00
102604	60 ml.	120 Pcs.	14.00	1650.00
102605	125 ml.	96 Pcs.	19.00	1728.00
102606	250 ml.	48 Pcs.	31.50	1464.00
102607	500 ml.	24 Pcs.	41.00	984.00
102608	1000 ml.	18 Pcs.	74.50	1350.00
102609	2000 ml.	06 Pcs.	200.00	1140.00
102610	4000 ml.	04 Pcs.	380.00	1360.00



14 Wide Mouth Reagent Bottle – HDPE

Material: High Density Poly Ethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent

Leak-proof & designed to protect light-sensitive liquids. Rigid with excellent chemical resistance to most acids, bases & alcohols. Ideal for storing, packaging & shipping. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
102703	30 ml.	120 Pcs.	11.50	1350.00
102704	60 ml.	120 Pcs.	13.75	1650.00
102705	125 ml.	96 Pcs.	19.00	1728.00
102706	250 ml.	48 Pcs.	31.50	1464.00
102707	500 ml.	24 Pcs.	40.50	984.00
102708	1000 ml.	18 Pcs.	74.00	1350.00
102709	2000 ml.	06 Pcs.	195.00	1140.00
102710	4000 ml.	04 Pcs.	370.00	1360.00



15 Wide Mouth Reagent Bottle – Amber

Material: High Density Poly Ethylene with PP Closure

HSN Code: 39233090

Colour: Amber, Opaque

Leak-proof & designed to protect light-sensitive liquids. Rigid with excellent chemical resistance to most acids, bases & alcohols. Ideal for storing, packaging & shipping. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
102803	30 ml.	120 Pcs.	13.50	1500.00
102804	60 ml.	120 Pcs.	15.50	1860.00
102805	125 ml.	96 Pcs.	22.00	1968.00
102806	250 ml.	48 Pcs.	36.50	1680.00
102807	500 ml.	24 Pcs.	46.50	1092.00
102808	1000 ml.	18 Pcs.	87.00	1548.00
102809	2000 ml.	06 Pcs.	215.00	1140.00
102810	4000 ml.	04 Pcs.	410.00	1360.00



16 Wide Mouth Reagent Bottle – LDPE

Material: Low Density Polyethylene with PP Closure

HSN Code: 39233090

Colour: Natural White, Translucent

These translucent bottles are leakproof having. Chemical Resistance, excellent flexibility yet can resist impact. Contents are easily visible. These Bottles are not autoclavable. Suitable for freezer storage down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
102903	30 ml.	120 Pcs.	13.50	1350.00
102904	60 ml.	120 Pcs.	15.50	1650.00
102905	125 ml.	96 Pcs.	21.00	1728.00
102906	250 ml.	48 Pcs.	34.50	1464.00
102907	500 ml.	24 Pcs.	44.00	984.00
102908	1000 ml.	18 Pcs.	77.00	1350.00
102909	2000 ml.	06 Pcs.	220.00	1170.00
102910	4000 ml.	04 Pcs.	425.00	1380.00





19 Sterile Wide Mouth Bottle – HDPE (Individually Wrapped)

Material: High Density Poly Ethylene with PP Closure
Colour: Natural White, Transparent

HSN Code: 39233090
MOQ: 4000 Nos of individual Size

Gamma sterilized and leakproof. It is rigid, offers excellent chemical resistance, and is suitable for freezer storage down to -100°C. Each container is individually wrapped in a single poly bag.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10263ST	30 ml.	900 Pcs.	26.00	23400.00
10264ST	60 ml.	502 Pcs.	28.00	14056.00
10265ST	125 ml.	264 Pcs.	34.00	8976.00
10266ST	250 ml.	144 Pcs.	50.00	7200.00
10267ST	500 ml.	78 Pcs.	68.00	5304.00
10268ST	1000 ml.	36 Pcs.	120.00	4320.00



20 Sterile Wide Mouth Bottle - Amber (Individually Wrapped)

Material: High Density Poly Ethylene with PP Closure
Colour: Amber, Opaque

HSN Code: 39233090
MOQ: 4000 Nos of individual Size

Gamma sterilized and leakproof, this container protects light-sensitive liquids. It is rigid, offers excellent chemical resistance, and is suitable for freezer storage down to -100°C. Each container is individually wrapped in a single poly bag.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10283STA	30 ml.	900 Pcs.	29.00	26100.00
10284STA	60 ml.	502 Pcs.	32.00	16064.00
10285STA	125 ml.	264 Pcs.	38.00	10032.00
10286STA	250 ml.	144 Pcs.	57.00	8208.00
10287STA	500 ml.	78 Pcs.	76.00	5928.00
10288STA	1000 ml.	36 Pcs.	132.00	4752.00



21 Wide Mouth Reagent Bottle – Print Graduation

Material: Polypropylene with PP Closure (In Conformance with US FDA 21 CF)
Colour: Natural, Transparent

HSN Code: 39233090

The round media storage wide-mouth bottles are made from durable material and feature clear printed scales for easy and precise observation, enhancing convenience. These bottles are leakproof, rigid, transparent, and offer excellent contact clarity. They are also autoclavable.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10381P	30 ml.	120 Pcs.	12.50	1500.00
10382P	60 ml.	120 Pcs.	15.00	1800.00
10383P	125 ml.	24 Pcs.	20.00	480.00
10384P	250 ml.	12 Pcs.	55.00	408.00
10385P	500 ml.	12 Pcs.	45.00	540.00
10386P	1000 ml.	06 Pcs.	75.00	450.00



22 Narrow Mouth Reagent Bottle – Print Graduation

Material: Polypropylene with PP Closure (In Conformance with US FDA 21 CF)
Colour: Natural, Transparent

HSN Code: 39233090

The round media storage Narrow-mouth bottles are made from durable material and feature clear printed scales for easy and precise observation, enhancing convenience. These bottles are leakproof, rigid, translucent, and are also autoclavable.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10351P	8 ml.	120 Pcs.	7.00	840.00
10352P	15 ml.	120 Pcs.	8.00	960.00
10353P	30 ml.	120 Pcs.	11.00	1320.00
10354P	60 ml.	120 Pcs.	14.00	1680.00
10355P	125 ml.	96 Pcs.	18.50	1776.00
10356P	250 ml.	48 Pcs.	29.50	1416.00
10357P	500 ml.	24 Pcs.	37.00	888.00
10358P	1000 ml.	18 Pcs.	71.00	1278.00

23 Conical Flask – Print Graduation

Material: Polypropylene with PP Closure
Colour: Natural, Transparent

HSN Code: 39233090

These transparent conical bottles are leakproof having excellent flexibility yet can resist impact. Contents are easily visible. These Bottles are autoclavable.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
103001	125 ml.	12 Pcs.	43.00	516.00
103002	250 ml.	12 Pcs.	62.00	744.00
103003	500 ml.	12 Pcs.	100.00	1200.00



24 Conical Flask – Graduation Engraved

Material: Polypropylene with PP Closure
Colour: Natural, Transparent

HSN Code: 39233090

Leakproof conical flask, also known as an Erlenmeyer flask has a broad, flat base and a conical body with a narrow neck are marked with clear, visible graduations for approximate volume measurements.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
103006	125 ml.	12 Pcs.	37.00	468.00
103007	250 ml.	12 Pcs.	51.00	588.00
103008	500 ml.	12 Pcs.	65.00	780.00



25 Volumetric Flask – Print Graduation

Material: Polypropylene with PP Closure
Colour: Natural, Transparent

HSN Code: 39233090

Class B volumetric flasks are widely used in laboratories for precise measurement of liquid volumes. Leakproof, Graduated with a single, precise graduation mark for accurate measurement.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
104101	25 ml.	06 Pcs.	52.00	312.00
104102	50 ml.	06 Pcs.	74.00	444.00
104103	100 ml.	06 Pcs.	91.00	546.00
104104	250 ml.	02 Pcs.	114.00	228.00
104105	500 ml.	02 Pcs.	162.00	324.00
104106	1000 ml.	02 Pcs.	240.00	480.00



26 Rectangular Wide Mouth Bottle - HDPE

Material: High Density Poly Ethylene with PP Closure
Colour: Natural White, Translucent

HSN Code: 39233090

Rigid and excellent chemical resistance to most acids, bases, and alcohols. Save space. Provides a Large labelling surface. Ideal for storing, packaging, and shipping dry materials or liquids. Good for freezer storage to down to -100°C.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10535H	125 ml.	12 Pcs.	20.00	240.00
10536H	250 ml.	12 Pcs.	34.00	408.00
10537H	500 ml.	12 Pcs.	44.00	528.00
10538H	1000 ml.	12 Pcs.	78.00	936.00
10539H	2000 ml.	06 Pcs.	225.00	1350.00





27 Rectangular Wide Mouth Bottle - Amber

Material: High Density Poly Ethylene with PP Closure
Colour: Amber, Opaque

HSN Code: 39233090

Rigid and excellent chemical resistance to most acids, bases, and alcohols. Save space. Provides a Large labelling surface. Ideal for storing, packaging, and shipping dry materials or liquids. Good for freezer storage to down to -100°C

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10535A	125 ml.	12 Pcs.	23.00	276.00
10536A	250 ml.	12 Pcs.	38.00	456.00
10537A	500 ml.	12 Pcs.	48.00	576.00
10538A	1000 ml.	12 Pcs.	90.00	1080.00
10539A	2000 ml.	06 Pcs.	240.00	1440.00



28 Large Wide Mouth Bottle with Handle – PP

Material: Polypropylene with PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Autoclavable. Ideal for storing and sterilizing contaminated lab items. Leakproof, space-saving design, wide mouth for easy access, integrated handle for better grip and pouring convenience, suitable for low-temperature applications (down to -100°C/-148°F).

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10371PR	4 Liters / 1 Gal. Round	06 Pcs.	450.00	2700.00
10372PS	4 Liters / 1 Gal. Square	06 Pcs.	450.00	2700.00



29 Large Wide Mouth Bottle with Handle – HDPE

Material: High Density Poly Ethylene with PP Closure
Colour: Natural White, Translucent

HSN Code: 39233090

Autoclavable. Ideal for storing and sterilizing contaminated lab items. Leakproof, space-saving design, wide mouth for easy access, integrated handle for better grip and pouring convenience, suitable for low-temperature applications (down to -100°C/-148°F).

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10373HR	4 Liters / 1 Gal. Round	06 Pcs.	450.00	2700.00
10374HS	4 Liters / 1 Gal. Square	06 Pcs.	450.00	2700.00



30 Large Narrow Mouth Bottle – PP

Material: Polypropylene with PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Leakproof & Autoclavable. Ideal for storing chemicals, standards, and distilled water. 4L and 5L have built-in shoulder for an identification tag.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
10375PP	02 Liters	06 Pcs.	230.00	1380.00
10376PP	04 Liters	01 Pcs.	420.00	420.00
10377PP	08 Liters	01 Pcs.	780.00	780.00

31 Aspirator Bottles with Stopcock, PP

Material: Polypropylene with PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Leakproof & Autoclavable. Ideal for storing and dispensing lab solutions and sterile / distilled water. The spigot itself can be autoclaved in the open position. Spigot is made of Styrene and Silicon. Note: Remove the spigot before autoclaving.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101114	02 Liters	01 Pcs.	620.00	620.00
101115	05 Liters	01 Pcs.	850.00	850.00
101116	10 Liters	01 Pcs.	1400.00	1400.00
101117	20 Liters	01 Pcs.	2000.00	2000.00
101119	Spare Spigot	01 Pcs.	300.00	300.00



32 Carboys, PP

Material: Polypropylene with white PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Graduated & Autoclavable, Leakproof assurance, shatterproof construction, excellent chemical resistance, suitable for storing, packaging, and shipping liquids, low in trace elements. Wide shoulder handle for easy pouring & carrying.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101121	05 Liters	01 Pcs.	800.00	850.00
101122	10 Liters	01 Pcs.	1200.00	1200.00
101123	20 Liters	01 Pcs.	1480.00	1480.00



33 Carboys Wide Mouth, PP

Material: Polypropylene with white PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Autoclavable & Graduated Wide Mouth for easy filling/dispensing for solid powder etc. Leakproof guarantee, shatter-resistant, suitable for bulk media and chemical storage.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101125	10 Liters	01 Pcs.	1350.00	1350.00
101126	20 Liters	01 Pcs.	2000.00	2000.00



34 Carboys With Stopcock, PP

Material: PP with white PP Closure (Confirming to US FDA 21CFR)
Colour: Natural, Transparent

HSN Code: 39233090

Autoclavable, Graduated, Rugged construction, compatible with most lab reagents including acids, bases, alcohols, aldehydes, esters, and aliphatic hydrocarbons, suitable for higher temperature applications up to 135°C. Preassembled leak-proof stopcock. Note: Remove the spigot before autoclaving.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
101131	10 Liters	01 Pcs.	1400.00	1400.00
101132	20 Liters	01 Pcs.	2000.00	2000.00
101133	Spare Spigot	01 Pcs.	300.00	300.00





34 Chemical Wide Mouth Bottles Amida shape - HDPE

Material: High Density Polyethylene with PP Closure
 Colour: Natural White, Translucent

HSN Code: 39233090

These Chemical bottles are Suitable for storing a variety of chemicals, including cleaners, detergents, oils, and liquid chemicals. Secure sealing options to prevent leakage

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
103601	50 ml.	120 Pcs.	13.00	1560.00
103602	100 mL.	96 Pcs.	19.00	1824.00
103603	250 mL.	48 Pcs.	32.00	1536.00
103604	500 mL.	24 Pcs.	40.00	960.00
103605	1000 mL.	18 Pcs.	74.00	1350.00



35 Oblong Jar with Handle (Rectangular Shaped)

Material: Polypropylene (P.P), Autoclavable
 Colour: Natural, Transparent

HSN Code: 39269099

These space-efficient, rectangular-shaped jars typically feature a wide mouth and are equipped with a handle for easy filling. They are ideal for transporting, pouring, and dispensing bulk items such as lab reagents, powders, and other dry goods. The large label panels make them perfect for adding product labels or branding.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
104008	2 Liters	06 Pcs.	190.00	1020.00



36 Squeeze (Gel) Bottle Dispenser

Material: Low Density Polyethylene (L.D.P.E.), Autoclavable
 Colour: Natural, Transparent

HSN Code: 39269099

A squeeze (gel) bottle dispenser is a versatile bottle commonly used in laboratories and various industries for precise dispensing of liquids. including reagents and buffers. Easy to clean and maintain, ensuring longevity and reliability.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
104001	225 mL. / 8 Oz.	24 Pcs.	25.50	612.00
104002	475 mL. / 16 Oz.	12 Pcs.	34.00	408.00
104003	750 mL. / 25 Oz.	12 Pcs.	42.00	504.00



37 Pharmaceutical Powder & Tablet Container (Multiuse)

Material: Polypropylene (P.P), Autoclavable

HSN Code: 39269099

Colour: Natural White, Translucent (Closures available in multi colour)

These bottles are suitable to store Gel, Tablets, Chemical Powder, Sterilized needles.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
103301	50 mL.	12 Pcs.	12.00	144.00
103302	100 mL.	12 Pcs.	17.00	204.00
103303	250 mL.	12 Pcs.	26.00	312.00

39 Measuring Cylinder (Hex Base)

Material: Polypropylene with PP Closure

HSN Code: 39269099

Colour: Natural, Transparent

These moulded graduations provide precise and easy-to-read measurements. Designed to withstand temperatures up to 60°C without permanent deformation, the stable base prevents tipping. Perfect for measuring and mixing liquids, solutions, and reagents in laboratory environments.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
113601H	10 ml.	24 Pcs.	13.00	312.00
113602H	25 ml.	12 Pcs.	19.50	234.00
113603H	50 ml.	12 Pcs.	23.50	282.00
113604H	100 ml.	12 Pcs.	27.50	330.00
113605H	250 ml.	12 Pcs.	54.00	648.00
113606H	500 ml.	12 Pcs.	81.00	972.00
113607H	1000 ml.	06 Pcs.	180.00	1080.00



40 Measuring Cylinder (Round Base)

Material: Polypropylene with PP Closure

HSN Code: 39269099

Colour: Natural, Transparent

These clear, graduated Polypropylene Cylinders feature a round base for increased stability. They can endure temperatures up to 60°C without exceeding tolerance limits permanently. The fine graduations provide easy visibility, making them perfect for measuring, mixing solutions & reagents in laboratory.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
113602R	25 ml.	12 Pcs.	20.00	240.00
113603R	50 ml.	12 Pcs.	24.00	288.00
113605R	250 ml.	12 Pcs.	53.00	636.00
113606R	500 ml.	12 Pcs.	82.00	984.00
113607R	1000 ml.	06 Pcs.	150.00	900.00
113608R	2000 ml.	01 Pcs.	500.00	500.00



41 Measuring Jug

Material: Polypropylene with PP Closure

HSN Code: 39269099

Colour: Natural, Transparent

Optically clear with easy-to-read, raised graduations. Tapered spout for easy pouring and transferring of liquids. Handle for safe and convenient lifting. Resistant to strong acids, alcohols, and bases

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
113902	250 ml.	12 Pcs.	26.50	318.00
113903	500 ml.	12 Pcs.	35.00	420.00
113904	1000 ml.	06 Pcs.	53.50	321.00
113905	2000 ml.	06 Pcs.	94.50	567.00
113907	3000 ml.	03 Pcs.	175.00	525.00
113908	5000 ml.	01 Pcs.	270.00	270.00



42 Measuring Beaker

Material: High Density Poly Ethylene with PP Closure

HSN Code: 39269099

Colour: Natural White, Translucent

Rigid and excellent chemical resistance to most acids, bases, and alcohols. Save space. Provides a Large labelling surface. Ideal for storing, packaging, and shipping dry materials or liquids. Good for freezer storage to down to -100°C

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
1014201	50 ml.	24 Pcs.	8.50	204.00
1014202	100 ml.	12 Pcs.	9.50	114.00
1014203	250 ml.	12 Pcs.	19.00	228.00
1014204	600 ml.	12 Pcs.	31.00	372.00
1014205	1000 ml.	06 Pcs.	49.00	294.00
1014206	2000 ml.	06 Pcs.	100.00	600.00
1014207	5000 ml.	03 Pcs.	270.00	780.00





43

43 Measuring Jugs - Printed Graduation

Material: Polypropylene (P.P), Autoclavable
Colour: Natural, Transparent

HSN Code: 39269099

These Measuring Jugs have easy to read silk-screened graduations and are Translucent Autoclavable. The handle type is easy to use for safe handling of acids, alkalise, and other reagent chemical solutions.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
114002	250 ml.	12 Pcs.	38.00	456.00
114003	500 ml.	12 Pcs.	60.00	720.00
114004	1000 ml.	06 Pcs.	92.00	552.00



44

44 Measuring Beaker – Printed Graduation

Material: Polypropylene (P.P), Autoclavable
Colour: Natural, Transparent

HSN Code: 39269099

These Beakers have Easy to read silk-screened graduations. Translucent, Autoclavable and stackable. Have excellent clarity, strength and chemical resistance.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
114301	50 ml.	24 Pcs.	15.00	360.00
114302	100 ml.	24 Pcs.	16.00	284.00
114303	250 ml.	12 Pcs.	27.00	324.00
114304	600 ml.	06 Pcs.	45.00	270.00
114305	1000 ml.	06 Pcs.	82.00	492.00
114306	2000 ml.	06 Pcs.	140.00	840.00



45

45 Rain Gauze Cylinder

Material: Polypropylene (Transparent) or Polystyrene (Clear)
Colour: Natural

HSN Code: 39269099

These Rain Gauze Cylinders are available in clear transparent (Polystyrene) and clear translucent (Polypropylene) varieties, with a size of 200 cm², perfect for gauging rainwater

Code	Capacity	Material	Packing	Price / Piece (₹)	Price/ Pack (₹)
114501	200 ml.	P.P.	24 Pcs.	47.00	1128.00
114502	200 ml.	P.S.	24 Pcs.	55.00	1320.00



46

46 Rain Gauge Collector with Cylinder

Material: Polystyrene (P.S) and Polypropylene (P.P.)
Colour: Natural, Transparent

HSN Code: 39269099

These translucent outer body ensure a leak-proof tipping container featuring a wide mouth for easy collection of Rain Water. A cylinder moulded with Polystyrene is placed inside the container to measure the Rain Water.

Code	Capacity	Material	Packing	Price / Piece (₹)	Price/ Pack (₹)
114503	200 ml.	P.S./P.P.	01 Pcs.	250.00	250.00

47 Retort Stand (Burette Stand)

Material: Polystyrene & Polypropylene
Colour: White

HSN Code: 39269099

Non-Corrosive Retort Stand are designed to hold various clamps & rings that can secure flasks, test tubes, burettes & other laboratory wares. The vertical Rod measures 12.5 mm x 750 mm, attached to a sturdy, heavy base ensuring stability.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
125401	225 x 150 mm Side Hole	01 Pcs.	455.00	455.00
125402	300 x 225 mm Center Hole	01 Pcs.	780.00	780.00
125403	300 x 225 mm Side Hole	01 Pcs.	800.00	800.00



48 Burette

Material: Polymethyl Methacrylate / PMP
Colour: Natural, Transparent

HSN Code: 39269099

These burettes are designed for accurate and repeatable measurements, ideal for various laboratory applications. Each Burette is individually calibrated. Unaffected by dilute mineral acids and bases except for HF and NH₄OH. Not designed for alcohol and organic solvents. excellent for general use, schools, and industrial labs.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
125801	25 ml.	10 Pcs.	252.00	2520.00
125802	50 ml.	10 Pcs.	275.00	2750.00
125803	100 ml.	10 Pcs.	500.00	5000.00




49 Burette Clamp (Single/ Double)

Material: Polypropylene
Colour: White

HSN Code: 39269099

Non-corroding clamps hold burettes firmly in place. Tapered Wedges can firmly hold the Rods having diameter 1/2" or 3/8".

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124501	Single	12 Pcs.	52.00	624.00
124502	Double	12 Pcs.	104.00	1248.00



50 Burette Holder (with adjustable screw bolt)

Material: Polypropylene
Colour: White

HSN Code: 39269099

Designed for lateral mounting for excellent readability. Works well with varied diameters Rods and Burettes as well.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124503	With adjustable Screw Bolt	06 Pcs.	84.00	504.00





51 Funnel Holder

Material: Polypropylene
Colour: White

HSN Code: 39269099

Non-corroding funnel holder Non-Corrosive Funnel Holders with diameter 25 MM to 200 MM

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
124504	Single	12 Pcs.	58.00	696.00



52 Separatory Funnel Holder

Material: Polypropylene
Colour: White

HSN Code: 39269099

Non-corroding funnel holder with a front opening which provides a clear view of the solution and allows unobstructed placement of separating funnel

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
124507	Single	12 Pcs.	66.00	792.00




53 Separatory Funnel

Material: Polypropylene (P.P), Autoclavable
Colour: Natural, Transparent

HSN Code: 39269099

Excellent contact clarity and strength, Leakproof threaded screw cap, PP stopcock with a self-lubricating plug. Autoclavable for sterilization, meets laboratory standards.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
145101	100 ml.	02 Pcs.	200.00	400.00
145102	250 ml.	02 Pcs.	230.00	460.00
145103	500 ml.	02 Pcs.	270.00	540.00



54 Analytical Funnel

Material: Polypropylene (P.P), Autoclavable
Colour: Natural, Transparent

HSN Code: 39269099

Autoclavable. Excellent strength and chemical resistance. Designed for use with standard filter paper. The funnel body is at a 60° angle with internal ribs at 58° for rapid filtration.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
145301	35 mm.	108 Pcs.	5.00	540.00
145302	50 mm.	108 Pcs.	5.50	594.00
145303	62 mm.	108 Pcs.	10.00	1080.00
145304	75 mm.	72 Pcs.	11.00	792.00
145305	100 mm.	72 Pcs.	20.00	1440.00

55 Desiccator (Vacuum)

Material: Polypropylene (P.P), Polycarbonate (P.C.)
Colour: White (Base), Clear (Top)

HSN Code: 39269099

These thick-walled vacuum desiccators are crafted from high-quality plastics, making them lightweight & easy to handle. They can maintain a vacuum of 29 inches of mercury for 24 hours without greasing. The three-way PP stopcock features a Teflon key, and its port accommodates a 6.35 mm inner diameter tube for vacuum suction. Ideal for preserving moisture-sensitive products and achieving low humidity levels.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
205801	150 mm.	01 Pcs.	1365.00	1365.00
205802	200 mm.	01 Pcs.	2145.00	2145.00
205803	250 mm.	01 Pcs.	2755.00	2755.00
205804	300 mm.	01 Pcs.	4710.00	4710.00
2057SR	Silicon 'O'Ring for 150/200 mm.	01 Pcs.	350.00	350.00
2058SR	Silicon 'O'Ring for 250/300 mm.	01 Pcs.	400.00	400.00
2060SC	Spare Stopcock	01 Pcs.	620.00	620.00



56 Desiccator (Plain)

Material: Polypropylene (P.P), Polycarbonate (P.C.)
Colour: White (Base), Clear (Top)

HSN Code: 39269099

These desiccators are durable and efficient. Thick perforated polypropylene disc with a central hole for easy lifting. Ideal for preserving moisture-sensitive products. Not for vacuum applications.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
205601	150 mm.	01 Pcs.	900.00	900.00
205602	200 mm.	01 Pcs.	1500.00	1500.00
205603	250 mm.	01 Pcs.	2100.00	2100.00
205604	300 mm.	01 Pcs.	3400.00	3400.00
2057SR	Silicon 'O'Ring for 150/200 mm.	01 Pcs.	450.00	450.00
2058SR	Silicon 'O'Ring for 250/300 mm.	01 Pcs.	400.00	400.00
2060SC	Spare Stopcock	01 Pcs.	620.00	620.00



57 Petri Dish Autoclavable

Material: Polypropylene (P.P) (In Conformance with US FDA 21 CF)
Colour: Natural, Transparent

HSN Code: 39269099

Uniplast Petri Dish are Clear, Durable, reusable, autoclavable at 121°C, top-rack dishwasher safe. Typically, these comes non-sterile,

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
167101	50 mm.	36 Pcs.	07.50	120.00
167102	75 mm.	36 Pcs.	10.00	360.00
167103	100 mm.	36 Pcs.	12.00	468.00



58 Petri Dish Disposable (Sterilized)

Material: Polystyrene (P.S)
Colour: Natural, Clear

HSN Code: 39269099

MOQ for ETO or GAMMA: 20000 Nos

Uniplast Disposable Petri Dish offers Optimum flatness, high optical clarity, stacking ring and slip-ring features, sterilized by irradiation or EO, packaged in linear tear, medical style bags.

Dry Heat Sterile (Pathogen Free)

Code	Capacity	Packing Method	Packing	Price / Piece (₹)	Price/ Pack (₹)
167105P	90 x15 mm	Individually Packed	480 Pcs.	9.00	4320.00
167106P	90 x15 mm	Sleeve of 10 Pcs.	500 Pcs.	8.50	4250.00

ETO Sterile

Code	Capacity	Packing Method	Packing	Price / Piece (₹)	Price/ Pack (₹)
167107E	90 x15 mm	Individually Packed	480 Pcs.	10.50	5040.00
167108E	90 x15 mm	Sleeve of 10 Pcs.	500 Pcs.	10.00	5000.00

GAMMA Sterile

Code	Capacity	Packing Method	Packing	Price / Piece (₹)	Price/ Pack (₹)
167109G	90 x15 mm	Individually Packed	480 Pcs.	11.75	5640.00
167110G	90 x15 mm	Sleeve of 10 Pcs.	500 Pcs.	11.50	5750.00





59 Rack for Petri Dish

Material: Polymethyl Methacrylate (PMMA)
Colour: Natural, Transparent

HSN Code: 39269099

Threaded screw top for secure closure. Four grooves inside the jar hold 10 slides in total. Each groove holds two slides of the size 3" x 1" back-to-back.

Code	Capacity	Material	Packing	Price / Piece (₹)	Price / Pack (₹)
167111	60	For 90 mm Petri Dish	2 Pcs.	1400.00	2800.00
167112	56	For 60 mm Petri Dish	2 Pcs.	1275.00	2550.00




60 L Shaped Spreader

Material: Polystyrene (P.S.)
Colour: Available in multi colour

HSN Code: 39269099

The sterile smooth design allows for smooth and even spreading of samples across surfaces, such as agar plates, without gouging or cutting the medium. Commonly used in laboratories for spreading bacterial cultures or other liquid samples evenly on petri dishes or plates.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
227811	L-Spreader	100 Pcs.	7.75	7750.00




61 Soft Loop

Material: Polystyrene (P.S.)
Colour: Available in multi colour

HSN Code: 39269099

This non-toxic, disposable loop meets semi-quantitative standards for sample handling, with a specially treated surface to improve droplet adhesion.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
227808	Soft Loop - 1 µl	100 Pcs.	3.75	375.00
227809	Soft Loop - 10 µl	100 Pcs.	3.75	375.00
227810	Inoculation Needle	100 Pcs.	3.75	375.00



62 Stirring Rods

Material: Polypropylene (P.P.)
Colour: White

HSN Code: 39269099

Durable, lightweight, and chemical-resistant

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
227803	∅ 6 mm x H - 150 mm	12 Pcs.	6.00	72.00
227804	∅ 7 mm x H - 200 mm	12 Pcs.	9.50	114.00
227807	∅ 10 mm x H - 250 mm	12 Pcs.	12.50	150.00

63 Spatula

Material: Polymethyl Methacrylate (PMMA)

HSN Code: 39269099

Colour: Available in multi colour

Autoclavable. Double-ended design with a flat spatula on one end and a spoon on the other, suitable for measuring, spreading, mixing, and transferring chemicals. Durable, easy to clean, comfortable grip with Flat blade, ergonomic handle, heat-resistant. Excellent strength and chemical resistance.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
227802	150 ML	50 Pcs.	7.00	350.00
227803	200 ML	50 Pcs.	9.00	450.00



63

64 Measuring Scoop

Material: Polystyrene (P.S.)

HSN Code: 39269099

Colour: Available in multi colour

These autoclavable plastic measuring scoops are designed with Flat bottom for stability, stable handle, heat and chemical resistant. Perfect for sample collection, measuring and weighing in laboratories.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
227822	2 ml.	12 Pcs.	6.00	72.00
227822	5 ml.	12 Pcs.	8.00	96.00
227822	10 ml.	12 Pcs.	11.50	138.00
227822	25 ml.	12 Pcs.	15.00	180.00
227822	50 ml.	12 Pcs.	19.00	228.00
227822	100 ml.	12 Pcs.	27.00	324.00
227822	250 ml.	12 Pcs.	55.50	666.00
227822	500 ml.	06 Pcs.	84.00	504.00
227822	1000 ml.	06 Pcs.	99.00	599.00
227822	Set All Sizes	01 Pcs. Each Size	320.00	320.00



63

65 Slide Box (Detachable)

Material: Polystyrene (P.S.)

HSN Code: 39231090

Colour: Black or Grey

These Slide boxes made from Polystyrene are designed with removable lid. They are built to be durable, compact, and resistant to humidity. Indexed slots/grooves to keep 1" x 3" sized slides securely in place, and come with a numbered index for easy identification.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
167206	25 Slides	24 Pcs.	70.00	1680.00
167207	50 Slides	12 Pcs.	144.00	1728.00
167208	100 Slides	06 Pcs.	250.00	1500.00



64

66 Slide Box and Slide Mailer (Hinged)

Material: Polypropylene (P.P.)

HSN Code: 39231090

Colour: White

These Slide Boxes and Slide Mailer are designed with hinged lid and are suitable to hold 1" x 3" slides. The heavy walls of the boxes ensure that they do not warp, splinter, or crack. The base of the box has numbered index for easy identification. Widely used in laboratories, pathology departments and educational institutions.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
167210	5 Slides (Slide Mailer)	48 Pcs.	13.00	624.00
167211	25 Slides	24 Pcs.	27.00	648.00
167212	50 Slides	12 Pcs.	55.00	660.00
167213	100 Slides	12 Pcs.	133.00	1596.00



65



67 Coplin Jar

Material: Polypropylene (P.P.)
Colour: White

HSN Code: 39231090

Threaded screw top for secure closure. Four grooves inside the jar hold 10 slides in total. Each groove holds two slides of the size 3" x 1" back-to-back.

Code	Capacity	Packing	Price / Piece (₹)	Price / Pack (₹)
167201	10 Slides	12 Pcs.	31.00	372.00



68 Kidney Tray

Material: Polypropylene (P.P.)

HSN Code: 90181100

Unioplast Kidney-shaped tray is designed with a flattened profile for easy carrying, easy cleaning and Autoclaving, making them essential in any operating room.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
238411	150 mm / 06 inch	100 Pcs.	20.00	2000.00
238412	200 mm / 08 inch	100 Pcs.	25.00	2500.00
238413	250 mm / 10 inch	100 Pcs.	34.00	3400.00
238414	300 mm / 12 inch	100 Pcs.	66.00	6600.00



69 Bed Pan

Material: Polypropylene with PP Closure

HSN Code: 90181100

Plastic bed pans with lids and handles. Lower profile for easy use, wider seat for comfort, superior splash control.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
238514	Basic	50 Pcs.	110.00	5500.00
238515	Premium	05 Pcs.	310.00	1550.00
238516	Superior	02 Pcs.	375.00	750.00



70 Urine Pot (Male and Female)

Material: Polypropylene (P.P) or High-Density Polyethylene (HDPE)
Colour: Natural, Transparent

HSN Code: 39233090

Male and Female Polypropylene Autoclavable Urine Pot 800 ml. Another option of White colour pot is also available made of HDPE. Easy fit removable lid for hygiene and odourless use, easy carry handle, can be stored in standing or lying position, portable design, graduated markings individually packed.

Code	Description	Material	Packing	Price / Piece (₹)	Price/ Pack (₹)
10341MP	Male Shape	PP	80 Pcs.	40.00	3200.00
10342FP	Female Shape	PP	80 Pcs.	40.00	3200.00
10341MH	Male Shape	HDPE	80 Pcs.	38.00	3040.00
10342FH	Female Shape	HDPE	80 Pcs.	38.00	3040.00

71 Sputum Mug with Lid

Material: Polypropylene (P.P)

HSN Code: 90181100

Colour: Natural, Transparent

Uniplast Sputum Mug has Secure-fitting lids to prevent spills and contamination. Lightweight, cost-effective, and commonly used in hospitals.

Code	Capacity	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
238310	300 ml	Bulk Pack	50 Pcs.	60.00	3000.00



71

72 Sputum Specimen Container

Material: Polypropylene (P.P)

HSN Code: 90181100

Colour: Natural, Transparent

MOQ: 2000 Nos for individual Pack

Uniplast Sputum Specimen sample container offers leak-resistant design with screw-on cap for secure fit, plastic spoon for sample collection, writing label for sample identification. Sputum Box

Code	Capacity	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
238315	Sputum - 30 ml	Bulk Pack	100 Pcs.	13.00	600.00
238316	Sputum - 30 ml	Individual Wrap	100 Pcs.	15.50	635.00



72

73 Stool Specimen Container

Material: Polypropylene (P.P)

HSN Code: 90181100

Colour: Natural, Transparent

MOQ: 2000 Nos for individual Pack

Available in two sizes, these specimen containers are Durable, chemical-resistant, and commonly used for sample storage. Transparent for easy sample visibility. Screw-on caps for secure sealing and contamination prevention.

Code	Capacity	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
238311	Stool - 30 ml	Bulk Pack	100 Pcs.	6.50	350.00
238312	Stool - 30 ml	Individual Wrap	100 Pcs.	9.00	380.00
238313	Stool - 50 ml	Bulk Pack	100 Pcs.	9.00	475.00
238314	Stool - 50 ml	Individual Wrap	100 Pcs.	11.50	510.00



73

74 Sample Container for Urine

Material: Polypropylene (P.P)

HSN Code: 90181100

Colour: Natural, Transparent

MOQ: 2000 Nos for individual Pack

Available in various sizes, these disposable containers are ideal for collecting urine. The leak-proof cap ensures safe sample transfer. These containers are offered either individually packed or in bulk packing.

Code	Capacity	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
238201	30 ml	Bulk Pack	100 Pcs.	3.50	350.00
238202	30 ml	Individual Wrap	100 Pcs.	3.80	380.00
238203	50 ml	Bulk Pack	100 Pcs.	5.00	500.00
238204	50 ml	Individual Wrap	100 Pcs.	5.25	525.00
238205	60 ml	Bulk Pack	100 Pcs.	6.00	600.00
238206	60 ml	Individual Wrap	100 Pcs.	6.25	625.00
238207	100 ml.	Bulk Pack	100 Pcs.	13.00	650.00
238208	100 ml.	Individual Wrap	100 Pcs.	15.50	775.00
238209	120 ml.	Bulk Pack	100 Pcs.	15.50	775.00
238210	120 ml.	Individual Wrap	100 Pcs.	18.00	900.00



74



75 Specimen Jar / Gas Jar

Material: Polystyrene (P.S.)
Colour: Clear

HSN Code: 39231090

These gas jars are designed for laboratory use with Dimensions: 50 mm (OD) x 200 mm (Height). Tight fit cap for specimen protection, Heavy foot & ground flange, used for preserving different kinds of specimens. Precision engineering, advanced manufacturing and calibration techniques.

Code	Particulars	Packing	Price / Piece (₹)	Price / Pack (₹)
167301	Cylindrical Shape	12 Pcs.	48.00	576.00



76 Wide Mouth Specimen Tube (Straight Wall)

Material: Polystyrene (P.S.) with LDPE Snap Cap / PP Screw Cap
Colour: Natural, Clear

HSN Code: 39269099

Uniplast Straight Wall Specimen Tubes are ideal for storage of specimens and sample blocks. Disposable transparent flat bottom specimen tubes with option of "Snap-On Polyethylene cap" and Polypropylene Cap. Widely used in medical, clinical, scientific research, and educational institutions

Code	Description	Capacity	Packing	Price / Piece (₹)	Price/Pack (₹)
167416	LDPE Snap Cap	7 Drams/25 ml 25 Ø x 75 H mm	100 Pcs.	21.00	2100.00
167417	LDPE Snap Cap	10 Drams/40 ml 37 Ø x 50 H mm	100 Pcs.	24.00	2400.00
167418	PP Screw Cap External Thread	16 Drams/60 ml 31Ø x100 H mm	100 Pcs.	36.00	3600.00



77 Laboratory Jars with Push Caps

Material: Polypropylene (P.P.) with Screw cap
Colour: White

HSN Code: 39231090

Translucent autoclavable container with excellent chemical resistance. Closures are grooved for easy stacking. Widely used in medical, clinical, scientific research, and educational institutions.

Code	Description	Ø x H mm	Packing	Price / Piece (₹)	Price/ Pack (₹)
167425	125 ml. / 4 Oz.	60 x 65 mm.	12 Pcs.	17.00	204.00
167426	250 ml. / 8 Oz.	60 x 100 mm.	12 Pcs.	29.00	348.00
167427	450 ml. / 15 Oz.	60 x 160 mm.	12 Pcs.	36.00	432.00



78 Swab Stick

Material: Polypropylene (P.P.) Tube with Wooden Stick
Colour: White

HSN Code: 39262021/ 39262010

Disposable Swab Sticks (Pack of 10) 6"inch Natural Eco-friendly Cotton Swab with Sterilized Sticks in Tubes. Cotton swabs designed with pure healthy cotton and durable cotton head.

Code	Capacity	Sterile	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167418	Swab Stick - 6"	Non-Sterile	100 Pcs.	7.00	700.00
167419	Swab Stick - 6"	Sterile	100 Pcs.	9.00	900.00
167420	Swab Stick -6" Without Tube	Non-Sterile	500 Pcs.	0.85	425.00

79 Specimen Sample Cup

Material: Polystyrene (P.S.)
Colour: Natural, Transparent

HSN Code: 39231090

Polystyrene (PS) sample cups are designed for Hitachi as well as Backman analysers for haematology and coagulation analysis of whole blood specimen, biochemical analysis of serum specimen.

Code	Capacity	Analyser	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167415	0.5 ml.	Backman Cup	1000 Pcs.	1.85	1850.00
167416	1.5 ml.	Backman Cup	1000 Pcs.	2.00	2000.00
167417	2 ml.	Hitachi Cup	500 Pcs.	1.10	550.00



80 Ria Vial

Material: Polystyrene with PP Closure
Colour: Natural, Amber

HSN Code: 90181100

RIA (Radioimmunoassay) vials made of Polystyrene are Cylindrical with round bottom. DNase, RNase, and pyrogen-free. Suitable for RIA, coagulation, and Bacteriology.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
167408	12 x 75 mm. Clear	200 Pcs.	1.10	220.00
167409	12 x 75 mm. (Light Amber)	200 Pcs	1.40	280.00
167410	Snap Caps	500 Pcs	1.00	420.00



81 Screw Top Micro Centrifuge tube with O-Ring Screw Cap

Material: Polypropylene
Colour: Natural, Transparent

HSN Code: 39231090

These tubes are designed to ensure sample integrity and prevent leaks, making them ideal for various laboratory applications. Screw top microtubes with O-ring screw caps for leak-proof sealing.

Code	Capacity	Packing	Price / Piece (₹)	Price/Pack (₹)
167424	1.5 ml.	500 Pcs.	5.25	2625.00



82 Micro Centrifuge Tube

Material: Polypropylene
Colour: Natural, Transparent

HSN Code: 39231090

Flip-top snap cap for secure closure. Moulded graduation marks for accurate measurement. Autoclavable and freezable. Conical bottom for efficient sample collection.

Code	Capacity	Packing	Price / Piece (₹)	Price/Pack (₹)
167419	0.1 ml.	1000 Pcs.	0.75	750.00
167420	0.5 ml.	1000 Pcs.	0.60	600.00
167421	1.5 ml.	500 Pcs.	0.70	350.00
167422	2.0 ml.	500 Pcs.	0.80	400.00





83

83 Storage Vials (self-standing)

Material: Polypropylene
Colour: Natural, Transparent

HSN Code: 39231090
MOQ: 5,000 Nos for Sterile Pack

Excellent chemical resistance, High strength, nonpyrogenic, nontoxic, and leakproof design. External graduations with marking area, and packed in Sealed bags.

Code	Capacity	Sterile	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167414	2.0 ml.	Non-Sterile	500 Pcs.	2.50	1250.00
167415	5.0 ml.	Non-Sterile	500 Pcs.	5.00	2500.00
167414S	2.0 ml.	Sterile	500 Pcs.	3.75	1875.00
167415S	5.0 ml.	Sterile	500 Pcs.	6.00	3000.00



84

84 Cryo Vial

Material: Polypropylene with PP Closure
Colour: Natural, Transparent

HSN Code: 39231090
MOQ: 5,000 Nos for Sterile Pack

Leak-proof design with moulded O-ring. Round bottom, self-standing options. RNase, DNase, and pyrogen-free. Suitable for storage at temperatures as low as -196°C. Ideal to use for storage of PCR Reagents, Enzymes and other diagnostics and biomedical samples.

Code	Capacity	Sterile	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167411	1.8 ml.	Non-Sterile	500 Pcs.	4.00	2000.00
167412	4.8 ml.	Non-Sterile	500 Pcs.	7.00	3500.00
167411S	1.8 ml.	Sterile	500 Pcs.	5.00	2500.00
167412S	4.8 ml.	Sterile	500 Pcs.	8.00	4000.00



85

85 Centrifugal Tubes (Conical and Self-standing Bottom)

Material: Polypropylene (P.P.)
Colour: Natural, Transparent

HSN Code: 39231090
MOQ: 40,000 Nos for Sterile Pack

Excellent chemical resistance, High strength, nonpyrogenic, nontoxic, and leakproof design. External graduations with marking area, and packed in Sealed bags.

Conical Bottom Centrifugal Tubes

Code	Capacity	Sterile	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167414	15 ml.	Non-Sterile	1600 Pcs.	8.00	12800.00
167415	50 ml.	Non-Sterile	500 Pcs.	9.50	4750.00
167414S	15 ml.	Sterile	1600 Pcs.	11.50	18400.00
167415S	50 ml.	Sterile	500 Pcs.	13.50	6750.00

Self-Standing Bottom Centrifugal Tubes

Code	Capacity	Sterile	Packing	Price / Pcs. (₹)	Price/ Pack (₹)
167415	50 ml.	Non-Sterile	500 Pcs.	10.50	5250.00
167415S	50 ml.	Sterile	500 Pcs.	14.25	7125.00



86

86 Conical Centrifuge Tube Rack

Material: Polypropylene (R.P.P.)
Colour: Orange or White

HSN Code: 39269099

Easy assembled Racks held together with joints. 30 wells to hold 15 ml tubes measuring 17 mm Dia and 20 wells to suitably hold 50 ml tubes measuring 30 mm Dia. Alphanumeric indexing aids for quick and easy tube identification.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
124008	15ml Tube with 17 mm Ø 50 ml Tubes with 30 mm Ø	02 Pcs.	215.00	430.00

87 Cryo Box

Material: Polypropylene

HSN Code: 39269099

Colour: Natural, White, Green

These Cryo boxes are designed to store vials at temperatures as low as -150°C. Autoclavable at +121°C for 20 minutes. Available in 10x10 (100 samples) configurations. Alpha-numeric coding, smart closing indication, clear lid, and coloured base.

Code	Particulars	Packing	Price / Piece (₹)	Price/ Pack (₹)
124210	100 Places for 1ml. or 1.8 ml.	04 Pcs.	225.00	900.00
124211	81 Places for 1ml. or 1.8 ml.	04 Pcs.	240.00	960.00



87

88 Multipurpose stand for MCT, RIA and Storage Vial

Material: Polypropylene (P.P.)

HSN Code: 39269099

Colour: White

These autoclavable racks feature 6 pegs designed to hold 12 mm diameter tubes (RIA vials, Cryo Vial) and 2 pegs suitable for 20 mm diameter Storage Vial. Convenient insertion and removal of tubes.

Code	Particulars	Packing	Price / Piece (₹)	Price/Packs (₹)
124003	Ø 12 mm & Ø 20 mm Tubes	12 Pcs.	45.00	540.00



88

89 Test Tube Stand (3 Tier)

Material: Polypropylene (R.P.P.)

HSN Code: 39269099

Colour: White

Test Tube Racks - Ingenious 3-tier design with convenient spaces for clear view of tube contents. Autoclavable and can withstand sub-freezing temperatures. Convenient insertion and removal of tubes. Available in various sizes and capacities.

Code	Places	Tube Ø Dia	Packing	Price / Piece (₹)	Price / Pack (₹)
124101	31 Places	13 mm	08 Pcs.	118.00	944.00
124103	31 Places	16 mm	08 Pcs.	118.00	944.00
124105	24 Places	18 mm	08 Pcs.	177.00	1416.00
124106	20 Places	20 mm	08 Pcs.	177.00	1416.00
124107	18 Places	25 mm	08 Pcs.	177.00	1416.00



89

90 Test Tube Stand


Material: Polypropylene (P.P.)

HSN Code: 39269099

Colour: White & Blue

These Test Tube stands moulded in polypropylene are designed to securely hold 6 test tubes (2 of 25 mm diameter and 4 of 16 mm diameter). Autoclavable, includes 6 vertical pins for drying tubes, and has corresponding hemispherical wells at the base.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124002	Ø 16 mm & Ø 25 mm	48 Pcs.	35.00	1680.00



90



91

91 Round Test Tube Stand

Material: Polypropylene (P.P.)
Colour: White

HSN Code: 39269099

These Polypropylene stands are designed to securely hold test tubes and are ideal for laboratory use. Holds 12 tubes (4 of 25mm diameter and 8 of 19 mm diameter). Space-saving, autoclavable, easy to disassemble for cleaning, includes 16 vertical pegs for drying test tubes, & a trough to contain drippings.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124001	Ø 19 mm & Ø 25 mm	12 Pcs.	154.00	1848.00



92

92 Pipette Stand (Vertical)

Material: Polypropylene (P.P.) Autoclavable
Colour: White

HSN Code: 39269099

Made of polypropylene, this stand is very stable, non-corrosive and detachable, holding 28 pipettes with a maximum diameter of 14 mm. It features a heavy base (225 /237 mm diameter) and an adjustable polypropylene-coated rod (375 /453 mm length). The top disc has holes for pipette placement, while the lower disc allows liquid drainage.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124201	Vertical Heavy Base	01 Pcs.	462.00	462.00
124202	Vertical Light Base (Economical)	06 Pcs.	335.00	2010.00



93


93 Pipette Rack (Horizontal)

Material: Polypropylene (P.P.) Autoclavable
Colour: White

HSN Code: 39269099

Made of **polypropylene**, this stand can hold **up to 12 pipettes**. It features three rods attaching the two sides firmly with strength. It is **autoclavable** and available in **white colour**.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124203	Horizontal	12 Pc.	145.00	1740.00



94

94 Draining Rack

Material: Acrylonitrile Butadiene Styrene/ Polycarbonate
Colour: White

HSN Code: 39269099

Features interchangeable ABS pegs, with closed-end sockets to prevent leakage. Includes a drip channel tray for waste water collection. 20 Peg Plastic Rack with 14 cm long pegs.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
124205	20 Pegs	01 Pc.	2150.00	2150.00
124206	42 Pegs	01 Pc.	4200.00	4200.00

95 Pipettor Stand (Z Type)

Material: Polymethyl Methacrylate (PMMA)

HSN Code: 39269099

A self-standing holder that accommodates most major brands of pipettes. It provides a stable and organized storage solution.

Code	Particulars	Packing	Price / Piece (₹)	Price/ Pack (₹)
205271	5 Places	02 Pcs.	440.00	880.00



95

96 Micro Pipette Stand with Drawer & Test Tube Rack

Material: Polymethyl Methacrylate (PMMA)

HSN Code: 39269099

Uniplast Micropipette Stand with Drawer is available in multiple variants to hold up to 3 or 6 Micropipettes. Can accommodate 22 tubes of 12/10 MM diameter. Comes with option of single and two separate drawers for keeping microtips.

Code	Particulars	Packing	Price / Piece (₹)	Price/ Pack (₹)
205273	3 PPT's Holder & 1 Drawer	02 Pcs.	525.00	1050.00
205274	6 PPT's Holder & 2 Drawers	02 Pcs.	775.00	1550.00
205275	3 PPT's Holder, 1 Drawers & 10 Places Test Tube Rack	02 Pcs.	575.00	1150.00
205276	6 PPT's Holder, 2 Drawers & 22 Places Test Tube Rack	02 Pcs.	850.00	1700.00



96

97 Micro Pipette – Single Channel – Fixed Volume

Material: PVDF / ABS, Autoclavable

HSN Code: 39269099

A lightweight, ergonomic micropipette, designed for fixed volume pipetting, this micropipette ensures high precision with minimum dead air space and soft plunger movement. It is fully autoclavable.

Code	Volume Range (µl)	Packing	Price / Piece (₹)	Price/ Pack (₹)
205209	2	01 Pcs.	1700.00	1700.00
205210	5	01 Pcs.	1700.00	1700.00
205211	10	01 Pcs.	1700.00	1700.00
205212	20	01 Pcs.	1700.00	1700.00
205213	25	01 Pcs.	1700.00	1700.00
205214	50	01 Pcs.	1700.00	1700.00
205215	100	01 Pcs.	1700.00	1700.00
205216	200	01 Pcs.	1700.00	1700.00
205217	250	01 Pcs.	1700.00	1700.00
205218	300	01 Pcs.	1700.00	1700.00
205219	400	01 Pcs.	1700.00	1700.00
205220	500	01 Pcs.	1700.00	1700.00
205221	1000	01 Pcs.	1700.00	1700.00



97

98 Micro Pipette – Single Channel – Variable Volume

Material: PVDF, ABS, Autoclavable

HSN Code: 39269099

A lightweight, ergonomic micropipette, designed for variable volume pipetting, this micropipette ensures high precision with minimum dead air space and soft plunger movement. It is fully autoclavable.

Code	Volume Range (µl)	Accuracy/Precision	Packing	Price/Pcs. (₹)	Price/Pack (₹)
205230	0.5 to 2	± 2.5% / <1.5%	01 Pcs.	2750.00	2750.00
205231	1 to 10	± 1.5% / <0.5%	01 Pcs.	2750.00	2750.00
205232	2 to 20	± 1.0% / <1.5%	01 Pcs.	2750.00	2750.00
205233	10 to 100	± 1.0% / <0.3%	01 Pcs.	2750.00	2750.00
205234	20 to 200	± 1.0% / <0.3%	01 Pcs.	2750.00	2750.00
205235	200 to 1000	± 1.0% / <0.2%	01 Pcs.	2750.00	2750.00
205236	1000 to 5000	± 0.5% / <0.2%	01 Pcs.	2750.00	2750.00
205237	1000 to 10000	± 0.5% / <0.1%	01 Pcs.	2750.00	2750.00



98



99 Bottle Top Dispenser

Material: PTFE (Piston, Valve, Discharge Assembly), Glass (Cylinder)

Colour: Grey

HSN Code: 39269099

Introducing the Uniplast Bottle Top Dispenser In accordance with ISO 8655 Standards., perfect for accurate and easy dispensing of most Acids, Bases and Solvents. Fully autoclavable at 121°C, 15 psi. Featuring a smooth plunger, adjustable delivery nozzle, leakproof Valve functions and springless PTFE valve manifold for high chemical compatibility. Glass Barrel is protected by PP Sleeve. Simple to use Calibration Functions. Easy to disassemble for cleaning and services.

Code	Volume Range (ml)	Increment	Accuracy	CV (±% ml)	Packing	Price/ Piece (Rs)
205240	0.25 to 2.5	0.05ml ±0.6%	±% 0.015 ml.	0.005	01 Pcs.	18000.00
205241	0.5 to 5	0.1 ml. ±0.5%	±% 0.025 ml.	0.005	01 Pcs.	18000.00
205242	1.0 to 10	0.2 ml. ±0.5%	±% 0.050 ml.	0.010	01 Pcs.	18000.00
205243	2.5 to 25	0.5 ml. ±0.5%	±% 0.125 ml.	0.025	01 Pcs.	21000.00
205244	2.5 to 30	0.5 ml. ±0.5%	±% 0.150 ml.	0.030	01 Pcs.	21000.00
205245	5.0 to 50	1.0 ml. ±0.5%	±% 0.250 ml.	0.050	01 Pcs.	21000.00
205246	5.0 to 60	1.0 ml. ±0.5%	±% 0.30 ml.	0.060	01 Pcs.	21000.00
205247	10 to 100	2.0 ml. ±0.5%	±% 0.50 ml.	0.100	01 Pcs.	26000.00



100 Micropipette Tips and Filter Tips

Material: Polypropylene (P.P.) Autoclavable

Colour: Natural, Blue, Yellow

HSN Code: 39269099

Code	Description (µl)	Packing	Price / Piece (₹)	Price/ Pack (₹)
205250	Tips - 10 µl (White)	1000 Pcs.	0.188	188.00
205251	Tips - 200 µl (White)	1000 Pcs.	0.184	184.00
205252	Tips - 200 µl (Yellow)	1000 Pcs.	0.235	235.00
205253	Tips - 1000 µl (White)	500 Pcs.	0.460	230.00
205254	Tips - 200 µl (Graduated)	1000 Pcs.	0.235	235.00
205255	Tips - 1000 µl (Graduated)	500 Pcs.	0.470	235.00
Filter Tips				
205256	Filter Tips - 10 µl	1000 Pcs.	1.790	1790.00
205257	Filter Tips - 20 µl	1000 Pcs.	1.910	1910.00
205258	Filter Tips - 100 µl	1000 Pcs.	2.025	2025.00
205259	Filter Tips - 200 µl	1000 Pcs.	2.140	2140.00
205260	Filter Tips -1000 µl	500 Pcs.	4.280	2140.00



101 Pipette Tip Box

Material: Polypropylene (P.P.) Autoclavable

HSN Code: 39231090

Uniplast Pipette Tip Boxes are manufactured from high-quality polypropylene, these boxes are autoclavable and transparent, providing a clear view of the tips inside.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
205260	96 Tips of 2 to 10 µl	06 Pcs.	160.00	960.00
205261	96 Tips of 2 to 200 µl	06 Pcs.	160.00	960.00
205262	100 Tips of 200 to 1000 µl	06 Pcs.	175.00	1050.00



102 Pipette Pump

Material: Acrylonitrile Butadiene Styrene/ Silicone/PP

HSN Code: 39269099

A fast-release pipetting device with a silicon rubber collet for holding various pipettes. It is acid- and alkali-resistant and can be disassembled for cleaning. Provides safe and accurate pipetting with a rotatable thumbwheel for liquid dispensing. Available in 2ml, 10ml, and 25ml variants.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
205201	2 ml.	01 Pc.	300.00	300.00
205202	10 ml.	01 Pc.	290.00	290.00
205203	25 ml.	01 Pc.	300.00	300.00


103 Pasture Pipette

Material: Low Density Poly Ethylene (L.D.P.E.)
Colour: Natural

HSN Code: 39269099

Graduated Disposable Pipettes, available from Uniplast. Made from LDPE material, these pipettes are Available in 1 ml., 3 ml, and 5 mL capacities, they are designed for precise liquid handling intended for single-use. Integral bulb design eliminates the need for plugging or autoclaving. Hydrophobic and inert, non-toxic to biological fluids and most acids.

Code	Particulars	Packing	Price / Piece (₹)	Price/ Pack (₹)
205204	1 ml.	500 Pcs.	1.60	800.00
205205	2 ml.	500 Pcs.	1.57	785.00
205206	5 ml.	500 Pcs.	2.90	1450.00



104 ESR Pipettes (Disposable)

Material: Optically clear Polystyrene

HSN Code: 39269099

Uniplast Easy-Rate disposable ESR system allows rapid execution of ESR tests using the Western Green method. These pipettes feature Fiber caps for automated liquid level adjustment and a printed graduated scale from 0 to 180.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
205208	ESR Pipettes Disposable – 3 ml	100 Pcs.	5.85	585.00



105 ESR Stand - For 12 Disposable ESR Pipette

Material: Acrylic

HSN Code: 39269099

A durable and convenient pipette holder with 12 precisely spaced holes designed to hold ESR pipettes exactly in a vertical position for accurate testing.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
205209	ESR Stand – 12 Pipettes	01 Pcs.	500.00	500.00



106 Laboratory Tray

Material: Polypropylene (P.P.)
Colour: Natural, Transparent

HSN Code: 39269099

These trays are Multipurpose, Steam Autoclavable, Tapered walls for nesting.

Code	Dimensions (mm.)	Capacity	Packing	Price/ Pcs. (₹)	Price/ Pack (₹)
228002	L 260 x W 210 x H 52	2.0 L	05 Pcs.	128.00	640.00
228003	L 320 x W 260 x H 72	4.5 L	05 Pcs.	150.00	750.00
228004	L 405 x W 300 x H 78	6.5 L	05 Pcs.	215.00	1075.00
228005	L 480 x W 350 x H 93	11.5 L	05 Pcs.	325.00	1625.00





107 Animal Cage

Material: Polypropylene / Stainless Steel
Colour: Translucent

HSN Code: 39269099

PP Animal Cage complete with SS Lid, built in food feeder along with PP Water Bottle fitted with SS Nozzle and silicon cork. Removable trays and easy-access doors for simple cleaning and maintenance. Metallic Top offers optimal ventilation and easy access to feed the species.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
103905	Rat Cage 422 x 290 x 180 mm	01 Pcs.	1950.00	1950.00
103906	Mice Cage 290 x 220 x 140 mm	01 Pcs.	890.00	890.00
103901	Feeder Water Bottle 125 ml.	06 Pcs.	137.00	822.00
103902	Feeder Water Bottle 250 ml.	06 Pcs.	156.00	936.00



108 Plantation Pots

Material: Polypropylene (P.P.)
Colour: White

HSN Code: 39269099

These Plantation Pots comes in three different size variants. Made of polypropylene with a downward extended rim for convenient lifting, even when filled. The bottom has holes for drainage. Drainage holes at the bottom to prevent waterlogging.

Code	Capacity	Packing	Price / Piece (₹)	Price/ Pack (₹)
227901	100 mm.	20 Pcs.	32.00	640.00
227902	150 mm.	20 Pcs.	60.00	1200.00
227903	200 mm.	10 Pcs.	110.00	1100.00



109 Kipps Apparatus

Material: High Density Polypropylene (H.D.P.E.)
Colour: White

HSN Code: 39269099

Unioplast Kipps apparatus is accurately designed and manufactured to avoid leakage while during the experiment to preparing and storing the hydrogen sulphide. This apparatus consists of three main portions—an upper chamber to hold the liquid reactant (e.g. acid), a middle chamber to contain the solid reactant (e.g. Iron Sulphide), and a lower chamber meant to collect the residual outlet after liquid reacts with solid. Available in 250 mL, 550 mL, and 1000 mL capacities.

Code	Particulars	Packing	Price / Piece (₹)	Price/Pack (₹)
156101	Kipps 250 ml.	01 Pcs.	880.00	880.00
156102	Kipps 500 ml.	01 Pcs.	1040.00	1040.00
156103	Kipps 1000 ml.	01 Pcs.	1180.00	1180.00
S15606	Spare Thistle Funnel Set	01 Pcs.	90.00	90.00
S15607	Spare Stop Cock	01 Pcs.	190.00	190.00
S15608	Spare Kipps Bottom Cap	01 Pcs.	40.00	40.00



110 Over Flow Jar

Material: Polypropylene (P.P.)
Colour: Natural White

HSN Code: 39269099

Unioplast Over Flow Vessel is designed to fit standard laboratory setups with an overflow spout, useful to measure the displacement of liquids, such as in Archimedes' principal experiments.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
113201	Cylindrical Shape	24 Pcs.	36.00	864.00


111 Paschal Law Apparatus

Material: Polypropylene (P.P.) / High Density Polyethylene (H.D.P.E.)
 Colour: Natural Clear with Blue Piston

HSN Code: 39269099

Uniplast Manufactured Paschal Law apparatus consists of a **spherical Ethylene bulb with multiple holes**, connected to a **50mL plastic syringe** via **flexible rubber tubing**. Paschal law is demonstrated When liquid is pushed through the syringe, it is ejected with equal force from multiple holes in the bulb.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
158907	Paschal Law Apparatus	06 Pcs.	250.00	1500.00



111


112 Pneumatic Trough

Material: Polypropylene (P.P.)
 Colour: Natural, Transparent

HSN Code: 39269099

These Trough, made of propylene, are autoclavable and can be used for a wide range of purposes including storage, sterilizing, and drying labware.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
228001	ø-200 mm x h-100 mm	04 Pcs.	113.00	452.00



112

113 Simple Cell Pot (Battery Jar)

Material: Polypropylene (P.P.)
 Colour: Natural, Transparent

HSN Code: 39269099

Simple Cell Pot is used for making simple cells and voltmeters etc. This pot has very good contact clarity and excellent chemical resistance. Being in Polypropylene, this pot is much chemical resistant and durable than the one used conventionally.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
228008	Simple Cell Pot	12 Pcs.	59.00	708.00



113

114 Leclanché Cell Pot

Material: Polypropylene (P.P.)
 Colour: Natural, Transparent

HSN Code: 39269099

This pot is made of Polypropylene, used in educational settings to demonstrate the principles of electrochemistry. Suitable for educational purposes and laboratory experiments

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
228009	Leclanché Cell Pot	12 Pcs.	70.00	840.00



114


115

115 Siphon Automatic Apparatus

Material: Nylon/ L.D.P.E. Bulb / Polypropylene (P.P.)
Colour: White

HSN Code: 39269099

This Siphon is leaking proof, comes with stop cock used to transfer/ control liquids. Acid and Alkaline Proof.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
158908	Siphon Apparatus	06 Pcs.	260.00	1560.00


116

116 T - Type Tubing Connectors

Material: Polypropylene (P.P.)
Colour: Natural

HSN Code: 39269099

Uniplast made Autoclavable T Tubing connectors have multiple barbs per arm providing additional contact between the fitting and tubing ensures a tight, secure connection between lengths of tubing having the same inner diameter (I.D.). Chemical-resistant polypropylene for use with a wide range of liquids or gases

Code	ID Size (mm.)	ID Size (inch)	Packing	Price/ Pcs. (₹)	Price/ Pack (₹)
227101	T - 6.4 mm.	1/4"	100 Pcs.	9.00	900.00
227102	T - 8.0 mm.	5/16"	100 Pcs.	11.00	1100.00
227103	T - 9.5 mm.	3/8"	100 Pcs.	15.00	1500.00


117

117 L - Type Tubing Connectors

Material: Polypropylene (P.P.)
Colour: Natural

HSN Code: 39269099

Polypropylene Resin. Autoclavable. Connectors have multiple barbs per arm providing additional contact between the fitting and the tubing for extra gripping power.

Code	ID Size (mm.)	ID Size (inch)	Packing	Price/ Pcs. (₹)	Price/ Pack (₹)
227106	L - 6.4 mm.	1/4"	100 Pcs.	9.00	900.00
227107	L - 8.0 mm.	5/16"	100 Pcs.	11.00	1100.00
227108	L - 9.5 mm.	3/8"	100 Pcs.	15.00	1500.00


118

118 Y - Type Tubing Connectors

Material: Polypropylene (P.P.)
Colour: Natural

HSN Code: 39269099

Uniplast Y Tubing connectors are Autoclavable. Ensures a tight, secure connection between lengths of tubing having the same inner diameter (I.D.) Maximize liquid flow with non-tapered hose barbed Y-Type Connectors, moulded of polypropylene resin.

Code	ID Size (mm.)	ID Size (inch)	Packing	Price/ Pcs. (₹)	Price/ Pack (₹)
227111	Y - 6.4 mm.	1/4"	100 Pcs.	9.00	900.00
227112	Y - 8.0 mm.	5/16"	100 Pcs.	11.00	1100.00
227113	Y - 9.5 mm.	3/8"	100 Pcs.	15.00	1500.00

119 Crystal Model Set


Material: Polypropylene

HSN Code: 90230010

Model kits for assembling various crystals lattices, which are in skeletal form. The lattices represent classic structures & have been built to a large enough size so that their different faces can be seen. Picture set is meant for Sodium Chloride.

Diamond Set	Sodium Chloride Set	Graphite Set
Black Jacks: 57	Red Jacks: 14	Red Jacks: 12
Connectors: 84	Black Jacks: 13	Black Jacks: 45
	Connectors: 54	Connectors: 69, Long Black Connectors: 8

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
158804	Diamond Set	01 Pcs.	185.00	185.00
158805	Sodium Chloride Set	01 Pcs.	185.00	185.00
158804	Graphite Set	01 Pcs.	185.00	185.00




120 Atomic Model Set

Material: Polypropylene

HSN Code: 90230010

Uniplast offers Atomic model in three variants Junior, Medium and Senior Set models. Junior set comes with 60 balls, Medium Set comes with 120 Balls and senior set comes with 360 Balls.

Code	Description	Packing	Price / Piece (₹)	Price/ Pack (₹)
158701	60 Balls Junior Set	01 Pcs.	250.00	250.00
158702	120 Balls Medium Set	01 Pcs.	350.00	350.00
158703	360 Balls Senior Set	01 Pcs.	750.00	750.00



CHEMICAL RESISTANCE FOR UNIPLAST PRODUCTS

CLASSES OF SUBSTANCES AT ROOM TEMPERATURE	HDPE	LDPE	PC	PP	PS
ACIDS, DILUTE OR WEAK	Green	Green	Green	Green	Green
ACIDS, STRONG AND CONCENTRATED	Green	Green	Orange	Green	Yellow
ALCOHOLS, ALIPHATIC	Green	Green	Blue	Green	Green
ALDEHYDES	Blue	Blue	Yellow	Blue	Orange
BASES	Green	Green	Orange	Green	Green
ESTERS	Blue	Blue	Orange	Blue	Orange
HYDROCARBONS, ALIPHATIC	Blue	Yellow	Yellow	Blue	Orange
HYDROCARBONS, AROMATIC	Blue	Yellow	Orange	Yellow	Orange
HYDROCARBONS, HALOGENATED	Yellow	Orange	Orange	Yellow	Orange
KETONES	Blue	Blue	Orange	Blue	Orange
OXIDIZING AGENTS, STRONG	Yellow	Yellow	Orange	Yellow	Orange

- 30 Days of constant exposure cause no damage. Plastics may even tolerate for years.
- Little or no damages after 30 days of constant exposure to the reagents.
- Some effect after 7 days of constant exposure to the reagent like crazing, cracking, loss of strength or discoloration.
- Not recommended for continuous use. Immediate damage may occur.

PHYSICAL PROPERTIES OF RAW MATERIAL USED FOR UNIPLAST PRODUCTS


Resin	Max use Temp°C	Brittleness Temp°C	Transparency							Flexibility	Water Absorption%
				Autoclaving	Gas	Dry Heat	Radiation	Disinfectant	Specific Gravity		
HDPE	120	-100	Translucent	No	Yes	No	Yes	Yes	0.95	Rigid	20.01
LDPE	80	-100	Translucent	No	Yes	No	Yes	Yes	0.92	Excel	<0.01
PC	135	-135	Clear	Yes+	Yes	No	No	Yes	1.20	Rigid	0.35
PP	135	0	Translucent	Yes	Yes	No	No	Yes	0.90	Rigid	<0.02
PS	90	20	Clear	No	Yes	No	Yes	Some	1.05	Rigid	0.05

Sterilization :

- Autoclaving (1210C, 15 psig for 20 minutes) Clean and rinse item with distilled water before autoclaving. Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperature.
- Always completely disengage threads before autoclaving.
- Gas - Ethylene oxide formaldehyde.
- Dry heat -1600C for 120 minutes.
- Disinfectant - Benzalkonium chloride, formalin, ethanol etc.
- Radiation - gamma irradiated at 2.5 Mard with unstabilized plastic. Sterilizing reduces mechanical strength, Do not use PC vessels for vaccum application if they have been autoclaved

TERMS & CONDITIONS OF SALE

- 1 This Catalog Includes all “Uniplast” Brand Products.
- 2 GST Will be charged extra as applicable.
3. The Price for items in bulk requirement will be Quoted on request.
4. All Prices are FOR NEW Delhi.
5. Uniplast Responsibility ceases once the cargo is delivered to the carrier/transporter.
6. Goods will be dispatched in standard packing or multiple thereof.
7. Capacity measures are for laboratory use only and not for commercial practice.
8. Uniplast reserves the right to change / modify any product without prior notice.
9. All disputes shall be settled with in the jurisdiction of New Delhi.
10. This Catalog is effective 1st April, 2025, & cancels all over previous Catalogs.



UNIPLAST

Marketing Office: BP 174, Shalimar Bagh (West), Delhi 110088, India

Production Facility: D 221, Sector 3, D.S.I.I.D.C., Bawana Ind. Area, New Delhi - 110039, India

Tele: +91-9971700758, +91-9971800784, +91-8700221715,
+91-7827065011, +91-11-45520870

Email: UNIPLAST@mail.com , UNIPLAST@uniplast.us

Website: Uniplast.org.in